

## MAP4K1 antibody

### Product Information

Catalog No.:	FNab04980
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Serine/threonine-protein kinase, which may play a role in the response to environmental stress. Appears to act upstream of the JUN N-terminal pathway. May play a role in hematopoietic lineage decisions and growth regulation. Able to autophosphorylate.

### Immunogen information

Immunogen:	mitogen-activated protein kinase kinase kinase kinase 1
Synonyms:	Mitogen-activated protein kinase kinase kinase kinase 1 Hematopoietic progenitor kinase MAPK/ERK kinase kinase kinase 1 (MEK kinase kinase 1, MEKKK 1) MAP4K1 HPK1
Observed MW:	90 kDa
Uniprot ID :	Q92918

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Raji cells were subjected to SDS PAGE followed by western blot with FNab04980(MAP4K1 antibody) at dilution of 1:1000