

MAPRE2 antibody

Product Information

Catalog No.:	FNab05002
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in microtubule polymerization, and spindle function by stabilizing microtubules and anchoring them at centrosomes. May play a role in cell migration(By similarity).

Immunogen information

Immunogen:	microtubule-associated protein, RP/EB family, member 2
Synonyms:	Microtubule-associated protein RP/EB family member 2 APC-binding protein EB2 End-binding protein 2 (EB2) MAPRE2 RP1
Observed MW:	39-42 kDa
Uniprot ID :	Q15555

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05002(MAPRE2 antibody) at dilution of 1:1000