

MARVELD3 antibody

Product Information

Catalog No.:	FNab05018
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


As a component of tight junctions, plays a role in paracellular ion conductivity.

Immunogen information

Immunogen:	MARVEL domain containing 3
Synonyms:	MARVEL domain-containing protein 3 MARVELD3 MRVLDC3
Observed MW:	40 kDa
Uniprot ID :	Q96A59

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab05018(MARVELD3 Antibody) at dilution of 1:1000