

MAT2A antibody

Product Information

Catalog No.:	FNab05027
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Catalyzes the formation of S-adenosylmethionine from methionine and ATP.

Immunogen information

Immunogen:	methionine adenosyltransferase II, alpha
Synonyms:	S-adenosylmethionine synthase isoform type-2 (AdoMet synthase 2) Methionine adenosyltransferase 2 (MAT 2) Methionine adenosyltransferase II (MAT-II) MAT2A AMS2 MATA2
Observed MW:	40-45 kDa
Uniprot ID :	P31153

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05027(MAT2A antibody) at dilution of 1:500