

MBTD1 antibody

Product Information

Catalog No.:	FNab05044
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Putative Polycomb group(PcG) protein. PcG proteins maintain the transcriptionally repressive state of genes, probably via a modification of chromatin, rendering it heritably changed in its expressibility(By similarity). Specifically binds to monomethylated and dimethylated 'Lys-20' on histone H4.

Immunogen information

Immunogen:	mbt domain containing 1
Synonyms:	MBT domain-containing protein 1 MBTD1
Observed MW:	65-71 kDa
Uniprot ID :	Q05BQ5

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05044(MBTD1 Antibody) at dilution of 1:600