

CSF1 antibody

Product Information

Catalog No.:	FNab05068
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	colony stimulating factor 1(macrophage)
Synonyms:	Macrophage colony-stimulating factor 1 (CSF-1, M-CSF, MCSF) Lanimostim Proteoglycan macrophage colony-stimulating factor (PG-M-CSF) Processed macrophage colony-stimulating factor 1 Macrophage colony-stimulating factor 1 43 kDa subunit CSF1
Observed MW:	70 kDa
Uniprot ID :	P09603

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



U-937 cells were subjected to SDS PAGE followed by western blot with FNab05068(MCSF antibody) at dilution of 1:200