

ME2 antibody

Product Information

Catalog No.:	FNab05080
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	malic enzyme 2, NAD(+)-dependent, mitochondrial
Synonyms:	NAD-dependent malic enzyme, mitochondrial (NAD-ME) Malic enzyme 2 ME2
Observed MW:	54 kDa, 65 kDa
Uniprot ID :	P23368

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab05080(ME2 antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05080(ME2 Antibody) at dilution of 1:600