

## MED20 antibody

### Product Information

Catalog No.:	FNab05092
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

Component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional preinitiation complex with RNA polymerase II and the general transcription factors.


### Immunogen information

Immunogen:	mediator complex subunit 20
Synonyms:	Mediator of RNA polymerase II transcription subunit 20 Mediator complex subunit 20 TRF-proximal protein homolog (hTRFP) MED20 TRFP
Observed MW:	23 kDa
Uniprot ID :	Q9H944

### Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200 - 1:2000

Image:



K562 cells were subjected to SDS PAGE followed by western blot with FNab05092(MED20 antibody) at dilution of 1:1000