

MAP2K4 antibody

Product Information

Catalog No.:	FNab05116
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	mitogen-activated protein kinase kinase 4
Synonyms:	Dual specificity mitogen-activated protein kinase kinase 4 (MAP kinase kinase 4, MAPKK 4) JNK-activating kinase 1 MAPK/ERK kinase 4 (MEK 4) SAPK/ERK kinase 1 (SEK1) Stress-activated protein kinase kinase 1 (SAPK kinase 1, SAPKK-1, SAPKK1) c-Jun N-terminal kinase kinase 1 (JNKK) MAP2K4 JNKK1 MEK4 MKK4 PRKMK4 SEK1 SERK1 SKK1
Observed MW:	44 kDa
Uniprot ID :	P45985

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab05116(MAP2K4 antibody) at dilution of 1:50



RAW264.7 cells were subjected to SDS PAGE followed by western blot with FNab05116(MAP2K4 antibody) at dilution of 1:600