

MAP3K2 antibody

Product Information

Catalog No.:	FNab05119
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	mitogen-activated protein kinase kinase kinase 2
Synonyms:	Mitogen-activated protein kinase kinase kinase 2 MAPK/ERK kinase kinase 2 (MEK kinase 2, MEKK 2) MAP3K2 MAPKKK2 MEKK2
Observed MW:	65 kDa
Uniprot ID :	Q9Y2U5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab05119(MEKK2 antibody) at dilution of 1:1000