

## MEPE antibody

### Product Information

Catalog No.:	FNab05129
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Promotes renal phosphate excretion and modulates mineralization.

### Immunogen information


Immunogen:	matrix extracellular phosphoglycoprotein
Synonyms:	Matrix extracellular phosphoglycoprotein Osteoblast/osteocyte factor 45 (OF45) Osteoregulin MEPE
Observed MW:	98 kDa
Uniprot ID :	Q9NQ76

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung slide using FNab05129(MEPE Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05129(MEPE antibody) at dilution of 1:1500