

NTMT1 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab05138 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information

| | |
|--------------|--|
| Immunogen: | methyltransferase like 11A |
| Synonyms: | N-terminal Xaa-Pro-Lys N-methyltransferase 1 Alpha N-terminal protein methyltransferase 1A Methyltransferase-like protein 11A N-terminal RCC1 methyltransferase X-Pro-Lys N-terminal protein methyltransferase 1A (NTM1A) N-terminal Xaa-Pro-Lys N-methyltransferase 1, N-terminally processed NTMT1 C9orf32 METTL11A NRMT NRMT1 |
| Observed MW: | 25 kDa |
| Uniprot ID : | Q9BV86 |

Application

| | |
|---------------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB |
| Recommended dilution: WB: | 1:500-1:2000 |
| Image: | |



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05138(METTL11A antibody) at dilution of 1:800