

METTL6 antibody

Product Information

Catalog No.:	FNab05142
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probable methyltransferase.

Immunogen information


Immunogen:	methyltransferase like 6
Synonyms:	tRNA N(3)-methylcytidine methyltransferase METTL6 Methyltransferase-like protein 6 (hMETTL6) METTL6
Observed MW:	30-32 kDa
Uniprot ID :	Q8TCB7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab05142(METTL6 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05142(METTL6 antibody) at dilution of 1:300