

MFAP4 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab05146 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


Could be involved in calcium-dependent cell adhesion or intercellular interactions. May contribute to the elastic fiber assembly and/or maintenance(PubMed:26601954).

Immunogen information


| | |
|--------------|---|
| Immunogen: | microfibrillar-associated protein 4 |
| Synonyms: | Microfibril-associated glycoprotein 4 MFAP4 |
| Observed MW: | 36 kDa |
| Uniprot ID : | P55083 |

Application

| | |
|-----------------------|-----------------------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:20-1:200 |
| Image: | |



Immunohistochemistry of paraffin-embedded human ovary using FNab05146(MFAP4 antibody) at dilution of 1:100



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab05146(MFAP4 antibody) at dilution of 1:600