

MDK antibody

Product Information

Catalog No.:	FNab05182
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	midkine(neurite growth-promoting factor 2)
Synonyms:	Midkine (MK) Amphiregulin-associated protein (ARAP) Midgestation and kidney protein Neurite outgrowth-promoting factor 2 Neurite outgrowth-promoting protein MDK MK1 NEGF2
Observed MW:	None
Uniprot ID :	P21741

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, IHC
Recommended dilution:	IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer slide using FNab05182(Midkine Antibody) at dilution of 1:50