

## MIER1 antibody

### Product Information

Catalog No.:	FNab05184
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Transcriptional repressor regulating the expression of a number of genes including SP1 target genes. Probably functions through recruitment of HDAC1 a histone deacetylase involved in chromatin silencing.

### Immunogen information


Immunogen:	mesoderm induction early response 1 homolog(Xenopus laevis)
Synonyms:	Mesoderm induction early response protein 1 (Early response 1, Er1, Mi-er1, hMi-er1) MIER1 KIAA1610
Observed MW:	82 kDa
Uniprot ID :	Q8N108

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using FNab05184(MIER1 Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab05184(MIER1 antibody) at dilution of 1:500