

MIER2 antibody

Product Information

Catalog No.:	FNab05185
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcriptional repressor.

Immunogen information

Immunogen:	mesoderm induction early response 1, family member 2
Synonyms:	Mesoderm induction early response protein 2 (Mier2) MIER2 KIAA1193
Observed MW:	48-50 kDa
Uniprot ID :	Q8N344

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05185(MIER2 antibody) at dilution of 1:500