

## MIPEP antibody

### Product Information

Catalog No.:	FNab05195
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Cleaves proteins, imported into the mitochondrion, to their mature size.

### Immunogen information


Immunogen:	mitochondrial intermediate peptidase
Synonyms:	Mitochondrial intermediate peptidase (MIP) MIPEP MIP
Observed MW:	70-80 kDa
Uniprot ID :	Q99797

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab05195(MIPEP antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05195(MIPEP antibody) at dilution of 1:1000