

IMMT antibody

Product Information

Catalog No.:	FNab05199
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	inner membrane protein, mitochondrial(mitofilin)
Synonyms:	MICOS complex subunit MIC60 Cell proliferation-inducing gene 4/52 protein Mitochondrial inner membrane protein Mitofilin p87/89 IMMT HMP MIC60 MINOS2
Observed MW:	87 kDa
Uniprot ID :	Q16891

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF, WB, IP, FC
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab05199 (IMMT antibody) at dilution of 1:100



IP Result of anti-Mitofilin (IP:FNab05199, 3ug; Detection:FNab05199 1:1000) with Raji cells lysate 3800ug.



hela cells were subjected to SDS PAGE followed by western blot with FNab05199 (IMMT antibody) at dilution of 1:1000