

## MAP3K10 antibody

### Product Information

Catalog No.:	FNab05214
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	mitogen-activated protein kinase kinase kinase 10
Synonyms:	Mitogen-activated protein kinase kinase kinase 10 Mixed lineage kinase 2 Protein kinase MST MAP3K10 MLK2 MST
Observed MW:	104-115 kDa
Uniprot ID :	Q02779

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using FNab05214 (MAP3K10 Antibody) at dilution of 1:50



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab05214 (MAP3K10 antibody) at dilution of 1:300