

## MAP3K11 antibody

### Product Information

Catalog No.:	FNab05215
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	mitogen-activated protein kinase kinase kinase 11
Synonyms:	Mitogen-activated protein kinase kinase kinase 11 Mixed lineage kinase 3 Src-homology 3 domain-containing proline-rich kinase MAP3K11 MLK3 PTK1 SPRK
Observed MW:	60 kDa
Uniprot ID :	Q16584

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab05215(MLK3 antibody) at dilution of 1:50



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab05215(MLK3 antibody) at dilution of 1:500