

MME antibody

Product Information

Catalog No.:	FNab05231
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	membrane metallo-endoropeptidase
Synonyms:	Neprilysin Atriopeptidase Common acute lymphocytic leukemia antigen (CALLA) Enkephalinase Neutral endopeptidase 24.11 (NEP, Neutral endopeptidase) Skin fibroblast elastase (SFE) MME EPN
Observed MW:	90-110 kDa
Uniprot ID :	P08473

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab05231(MME,CD10 antibody) at dilution of 1:1000