

MMP8 antibody

Product Information

Catalog No.:	FNab05246
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Can degrade fibrillar type I, II, and III collagens.

Immunogen information


Immunogen:	matrix metalloproteinase 8(neutrophil collagenase)
Synonyms:	Neutrophil collagenase Matrix metalloproteinase-8 (MMP-8) PMNL collagenase (PMNL-CL) MMP8 CLG1
Observed MW:	70 kDa
Uniprot ID :	P22894

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC : 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab05246(MMP8 antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05246(MMP8 antibody) at dilution of 1:500