

## UBE2V2 antibody

### Product Information

Catalog No.:	FNab05250
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	ubiquitin-conjugating enzyme E2 variant 2
Synonyms:	Ubiquitin-conjugating enzyme E2 variant 2 DDVit 1 Enterocyte differentiation-associated factor 1 (EDAF-1) Enterocyte differentiation-promoting factor 1 (EDPF-1) MMS2 homolog Vitamin D3-inducible protein UBE2V2 MMS2 UEV2
Observed MW:	16-20 kDa
Uniprot ID :	Q15819

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab05250(UBE2V2 antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab05250(UBE2V2 antibody) at dilution of 1:1500