

## MNS1 antibody

### Product Information

Catalog No.:	FNab05256
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in the control of meiotic division and germ cell differentiation through regulation of pairing and recombination during meiosis.

### Immunogen information


Immunogen:	meiosis-specific nuclear structural 1
Synonyms:	Meiosis-specific nuclear structural protein 1 MNS1
Observed MW:	55-60 kDa
Uniprot ID :	Q8NEH6

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab05256(MNS1 Antibody) at dilution of 1:50



rat testis tissue were subjected to SDS PAGE followed by western blot with FNab05256(MNS1 Antibody) at dilution of 1:1000