

MOB1A antibody

Product Information

Catalog No.:	FNab05258
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	MOB kinase activator 1A
Synonyms:	MOB kinase activator 1A Mob1 alpha (Mob1A) Mob1 homolog 1B Mps one binder kinase activator-like 1B MOB1A C2orf6 MOB4B MOBK1B MOBK1B
Observed MW:	25 kDa
Uniprot ID :	Q9H8S9

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab05258(MOBKL1A antibody) at dilution of 1:100



IP Result of anti-MOBKL1A (IP:FNab05258, 4ug; Detection:FNab05258 1:1000) with Jurkat cells lysate 2800ug.



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05258(MOBKL1A antibody) at dilution of 1:500