

## MOB1B antibody

### Product Information


Catalog No.:	FNab05259
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	MOB kinase activator 1B
Synonyms:	MOB kinase activator 1B Mob1 homolog 1A (Mob1A) Mob1B Mps one binder kinase activator-like 1A MOB4A MOBKL1A
Observed MW:	28 kDa
Uniprot ID :	Q7L9L4

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05259( MOB1B Antibody) at dilution of 1:600