

MOB4 antibody

Product Information

Catalog No.:	FNab05261
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	MOB1, Mps One Binder kinase activator-like 3
Synonyms:	MOB-like protein phocein 2C4D Class II mMOB1 Mob1 homolog 3 (Mob3) Mps one binder kinase activator-like 3 Preimplantation protein 3 MOB4 MOB3 MOBKL3 PHOCN PREI3
Observed MW:	24 kDa
Uniprot ID :	Q9Y3A3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab05261(MOBKL3 Antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05261(MOBKL3 antibody) at dilution of 1:200