

MORC1 antibody

Product Information

Catalog No.:	FNab05272
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for spermatogenesis.

Immunogen information

Immunogen:	MORC family CW-type zinc finger 1
Synonyms:	MORC family CW-type zinc finger protein 1 Cancer/testis antigen 33 (CT33) MORC1 MORC
Observed MW:	108 kDa
Uniprot ID :	Q86VD1

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab05272(MORC1 antibody) at dilution of 1:500