

MOSPD3 antibody

Product Information

Catalog No.:	FNab05276
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	motile sperm domain containing 3
Synonyms:	Motile sperm domain-containing protein 3 MOSPD3
Observed MW:	28 kDa
Uniprot ID :	O75425

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab05276(MOSPD3 antibody) at dilution of 1:600