

MPC2 antibody

Product Information

Catalog No.:	FNab05279
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Mediates the uptake of pyruvate into mitochondria.

Immunogen information


Immunogen:	brain protein 44
Synonyms:	Mitochondrial pyruvate carrier 2 Brain protein 44 MPC2 BRP44
Observed MW:	18-20 kDa
Uniprot ID :	O95563

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:50-1:500
Image:	



Immunofluorescent analysis of HepG2 cells using FNab05279 (MPC2 antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab05279(MPC2 antibody) at dilution of 1:1000