

## MPP1 antibody

### Product Information

Catalog No.:	FNab05287
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Essential regulator of neutrophil polarity. Regulates neutrophil polarization by regulating AKT1 phosphorylation through a mechanism that is independent of PIK3CG activity(By similarity).

### Immunogen information

Immunogen:	membrane protein, palmitoylated 1, 55kDa
Synonyms:	55 kDa erythrocyte membrane protein (p55) Membrane protein, palmitoylated 1 MPP1 DXS552E EMP55
Observed MW:	55 kDa
Uniprot ID :	Q00013

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05287(MPP1 antibody) at dilution of 1:1000