

MPV17L antibody

Product Information

Catalog No.:	FNab05298
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Isoform 1 participates in reactive oxygen species metabolism by up-or down-regulation of the genes of antioxidant enzymes.

Immunogen information

Immunogen:	MPV17 mitochondrial membrane protein-like
Synonyms:	Mpv17-like protein M-LP homolog (M-LPH) MPV17L
Observed MW:	10 kDa
Uniprot ID :	Q2QL34

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab05298(MPV17L antibody) at dilution of 1:1000