

MPZL2 antibody

Product Information

Catalog No.:	FNab05301
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Mediates homophilic cell-cell adhesion.

Immunogen information

Immunogen:	myelin protein zero-like 2
Synonyms:	Myelin protein zero-like protein 2 Epithelial V-like antigen 1 MPZL2 EVA EVA1
Observed MW:	24 kD and 32kd
Uniprot ID :	O60487

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05301(MPZL2 antibody) at dilution of 1:800