

ABCC2 antibody

Product Information

Catalog No.:	FNab05310
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ATP-binding cassette, sub-family C(CFTR/MRP), member 2
Synonyms:	ATP-binding cassette sub-family C member 2 Canalicular multidrug resistance protein Canalicular multispecific organic anion transporter 1 Multidrug resistance-associated protein 2 ABCC2 CMOAT CMOAT1 CMRP MRP2
Observed MW:	190-250 kDa
Uniprot ID :	Q92887

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using FNab05310(ABCC2 Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05310(ABCC2 Antibody) at dilution of 1:600