

## ABCC5 antibody

### Product Information

Catalog No.:	FNab05311
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	ATP-binding cassette, sub-family C(CFTR/MRP), member 5
Synonyms:	ATP-binding cassette sub-family C member 5 Multi-specific organic anion transporter C (MOAT-C) Multidrug resistance-associated protein 5 SMRP pABC11 ABCC5 MRP5
Observed MW:	None
Uniprot ID :	O15440

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab05311(MRP5 Antibody) at dilution of 1:50