

## MRPL10 antibody

### Product Information

Catalog No.:	FNab05313
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	mitochondrial ribosomal protein L10
Synonyms:	Large ribosomal subunit protein uL10m 39S ribosomal protein L10, mitochondrial (L10mt, MRP-L10) 39S ribosomal protein L8, mitochondrial (L8mt, MRP-L8) MRPL10 MRPL8 RPML8
Observed MW:	29 kDa
Uniprot ID :	Q7Z7H8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab05313(MRPL10 antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab05313(MRPL10 antibody) at dilution of 1:500