

MRPL19 antibody

Product Information

Catalog No.:	FNab05319
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	mitochondrial ribosomal protein L19
Synonyms:	Large ribosomal subunit protein bL19m 39S ribosomal protein L15, mitochondrial (L15mt, MRP-L15) 39S ribosomal protein L19, mitochondrial (L19mt, MRP-L19) MRPL19 KIAA0104 MRPL15
Observed MW:	34 kDa
Uniprot ID :	P49406

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab05319(MRPL19 antibody) at dilution of 1:100



human heart tissue were subjected to SDS PAGE followed by western blot with FNab05319(MRPL19 antibody) at dilution of 1:500