

MRPL27 antibody

Product Information

Catalog No.:	FNab05324
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	mitochondrial ribosomal protein L27
Synonyms:	Large ribosomal subunit protein bL27m 39S ribosomal protein L27, mitochondrial (L27mt, MRP-L27) MRPL27
Observed MW:	10 kDa
Uniprot ID :	Q9P0M9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skin cancer using FNab05324(MRPL27 antibody) at dilution of 1:50



A431 cells were subjected to SDS PAGE followed by western blot with FNab05324(MRPL27 antibody) at dilution of 1:300