

MRPL3 antibody

Product Information

Catalog No.:	FNab05326
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	mitochondrial ribosomal protein L3
Synonyms:	Large ribosomal subunit protein uL3m 39S ribosomal protein L3, mitochondrial (L3mt, MRP-L3) MRPL3 MRL3 RPML3
Observed MW:	39 kDa
Uniprot ID :	P09001

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab05326(MRPL3 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05326(MRPL3 antibody) at dilution of 1:1000