

## MRPL33 antibody

### Product Information

Catalog No.:	FNab05327
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	mitochondrial ribosomal protein L33
Synonyms:	Large ribosomal subunit protein bL33m 39S ribosomal protein L33, mitochondrial (L33mt, MRP-L33) MRPL33 C2orf1
Observed MW:	None
Uniprot ID :	O75394

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, IF
Recommended dilution:	IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab05327(MRPL33 antibody) at dilution of 1:50