

MRPL34 antibody

Product Information

Catalog No.:	FNab05328
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	mitochondrial ribosomal protein L34
Synonyms:	Large ribosomal subunit protein bL34m 39S ribosomal protein L34, mitochondrial (L34mt, MRP-L34) MRPL34
Observed MW:	10 kDa
Uniprot ID :	Q9BQ48

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer using FNab05328(MRPL34 antibody) at dilution of 1:100



human liver tissue were subjected to SDS PAGE followed by western blot with FNab05328(MRPL34 antibody) at dilution of 1:2000