

MRPL55 antibody

Product Information

Catalog No.:	FNab05340
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	mitochondrial ribosomal protein L55
Synonyms:	Large ribosomal subunit protein mL55 39S ribosomal protein L55, mitochondrial (L55mt, MRP-L55) Mitochondrial large ribosomal subunit protein bL31m MRPL55
Observed MW:	13 kDa
Uniprot ID :	Q7Z7F7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab05340(MRPL55 Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab05340(MRPL55 Antibody) at dilution of 1:600