

## MRPS12 antibody

### Product Information

Catalog No.:	FNab05345
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	mitochondrial ribosomal protein S12
Synonyms:	Small ribosomal subunit protein uS12m 28S ribosomal protein S12, mitochondrial (MRP-S12, S12mt) MT-RPS12 MRPS12 RPMS12 RPSM12
Observed MW:	15 kDa
Uniprot ID :	O15235

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab05345(MRPS12 antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05345(MRPS12 antibody) at dilution of 1:500