

MRPS28 antibody

Product Information

Catalog No.:	FNab05359
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	mitochondrial ribosomal protein S28
Synonyms:	Small ribosomal subunit protein bS1m 28S ribosomal protein S28, mitochondrial (MRP-S28, S28mt) 28S ribosomal protein S35, mitochondrial (MRP-S35, S35mt) MRPS28 MRPS35
Observed MW:	21 kDa
Uniprot ID :	Q9Y2Q9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab05359(MRPS28 antibody) at dilution of 1:400