

MRPS5 antibody

Product Information

Catalog No.:	FNab05364
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	mitochondrial ribosomal protein S5
Synonyms:	Small ribosomal subunit protein uS5m 28S ribosomal protein S5, mitochondrial (MRP-S5, S5mt) MRPS5
Observed MW:	42 kDa
Uniprot ID :	P82675

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05364(MRPS5 Antibody) at dilution of 1:600