

## MRRF antibody

### Product Information

Catalog No.:	FNab05366
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Responsible for the release of ribosomes from messenger RNA at the termination of protein biosynthesis. May increase the efficiency of translation by recycling ribosomes from one round of translation to another(By similarity).

### Immunogen information


Immunogen:	mitochondrial ribosome recycling factor
Synonyms:	Ribosome-recycling factor, mitochondrial (RRF, mtRRF) Ribosome-releasing factor, mitochondrial MRRF
Observed MW:	29 kDa
Uniprot ID :	Q96E11

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab05366(MRRF antibody) at dilution of 1:50



human lung tissue were subjected to SDS PAGE followed by western blot with FNab05366(MRRF antibody) at dilution of 1:500