

MSI2 antibody

Product Information

Catalog No.:	FNab05377
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


RNA binding protein that regulates the expression of target mRNAs at the translation level. May play a role in the proliferation and maintenance of stem cells in the central nervous system(By similarity).

Immunogen information


Immunogen:	musashi homolog 2
Synonyms:	RNA-binding protein Musashi homolog 2 (Musashi-2) MSI2
Observed MW:	37-40 kDa
Uniprot ID :	Q96DH6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab05377(Musashi-2 antibody) at dilution of 1:10



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab05377(Musashi-2 antibody) at dilution of 1:500