

MSRB2 antibody

Product Information


Catalog No.:	FNab05383
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	methionine sulfoxide reductase B2
Synonyms:	Methionine-R-sulfoxide reductase B2, mitochondrial (MsrB2) MSRB2 CBS-1 MSRB
Observed MW:	20 kDa
Uniprot ID :	Q9Y3D2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded human esophagus using FNa05383(MSRB2 Antibody) at dilution of 1:100



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNa05383(MSRB2 antibody) at dilution of 1:1000