

MTA3 antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab05390 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |


Background

Immunogen information


| | |
|--------------|--|
| Immunogen: | metastasis associated 1 family, member 3 |
| Synonyms: | Metastasis-associated protein MTA3 MTA3 KIAA1266 |
| Observed MW: | 65 kDa |
| Uniprot ID : | Q9BTC8 |

Application

| | |
|-----------------------|---------------------------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:500 - 1:2000; IHC: 1:50 - 1:200 |
| Image: | |



Immunohistochemistry of paraffin-embedded rat liver using FNab05390(MTA3 antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05390(MTA3 antibody) at dilution of 1:1000